



## Product description

|                     |                                                                                                                                                                                                                          |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 83.3942.501                                                                                                                                                                                                              |
| Product description | Filtropur V100, Vacuum filtration unit, 1,000 ml, screw cap, membrane: PES, Ø membraner: 91 mm, pore size: 0.22 µm, for sterile filtration, sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase, 1 piece(s)/bag |

## Product characteristics

|                        |                        |
|------------------------|------------------------|
| Graduation             | yes                    |
| Closure type           | screw cap              |
| Example of use         | for sterile filtration |
| Filter                 | yes                    |
| Surface area of filter | 58.93 cm <sup>2</sup>  |
| Pore size              | 0,22                   |
| Connection             | Step taper             |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### Size

|                   |             |
|-------------------|-------------|
| Sample volume     | 1.000.00 ml |
| Filtration volume | 1000 ml     |
| Height of product | 31 cm       |

### Material & colours

|                   |                                   |
|-------------------|-----------------------------------|
| Product material  | Polystyrene (PS)                  |
| Colour of product | transparent                       |
| Closure material  | High Density Polyethylene (HD-PE) |
| Membrane material | Polyethersulfone (PES)            |

### Purity & certification

|                  |                                                                 |
|------------------|-----------------------------------------------------------------|
| Product category | no medical device   IVD                                         |
| Purity standard  | sterile, non-pyrogenic, non-cytotoxic, free from<br>DNase/Rnase |
| Sterilisation    | Electron irradiation                                            |
| Non-cytotoxic    | yes                                                             |
| Batched          | yes                                                             |

### Packaging

|                                  |               |
|----------------------------------|---------------|
| Minimum order qty.               | 12            |
| Type of smallest<br>subpackaging | bag           |
| Piece(s) / inner box             | 1             |
| Piece (s) / outer case           | 12            |
| Piece(s) / pallet                | 192           |
| Depth of case                    | 545 mm        |
| Width of case                    | 425 mm        |
| Height of case                   | 350 mm        |
| Case volume                      | 0.0811 cbm    |
| Weight of product                | 0.26 kg       |
| Weight of case                   | 5.42 kg       |
| EAN of inner box                 | 4038917969123 |
| EAN case                         | 4038917969130 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

