



Product description

| | |
|---------------------|--|
| Order number | 83.3945.999 |
| Product description | Cell strainer adapter, for 15 ml, 5 ml, 17 x 100 mm and 12 x 75 mm tubes, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister |

Product characteristics

| | |
|--------------|---|
| Suitable for | 15 ml, 5 ml, 17 x 100 mm and 12 x 75 mm tubes |
|--------------|---|

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|--|
| Product category | no medical device IVD |
| Purity standard | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation | Gas sterilisation |
| Pyrogen-free | yes |
| Batched | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Cell strainer adapter, for 15 ml, 5 ml, 17 x 100 mm
and 12 x 75 mm tubes

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 25 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 25 |
| Piece (s) / outer case | 200 |
| Piece(s) / pallet | 8400 |
| Depth of case | 229 mm |
| Width of case | 330 mm |
| Height of case | 584 mm |
| Case volume | 0.0441 cbm |
| Weight of product | 0.0053 kg |
| Weight of case | 2 kg |
| EAN of inner box | 4038917415323 |
| EAN case | 4038917415316 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

