



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 65.386.992   |
| Product description | Colour-coded inserts, for CryoPure tubes, colour mix, 100 piece(s)/bag |

### **Product characteristics**

|              |                |
|--------------|----------------|
| Graduation   | no             |
| Suitable for | CryoPure tubes |

### **Material & colours**

|                   |                  |
|-------------------|------------------|
| Product material  | Polystyrene (PS) |
| Colour of product | colour mix       |

### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 2500          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 2500          |
| Piece(s) / pallet             | 810000        |
| Depth of case                 | 268 mm        |
| Width of case                 | 191 mm        |
| Height of case                | 100 mm        |
| Case volume                   | 0.0051 cbm    |
| Weight of product             | 0.0002 kg     |
| Weight of case                | 0.6 kg        |
| EAN of inner box              | 4038917228541 |
| EAN case                      | 4038917228534 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

