



**Product description**

|                     |  |
|---------------------|--|
| Order number        | 72.383.007   |
| Product description | CryoPure tubes, 5 ml, tube: PP, QuickSeal screw cap, cap assembled, HD-PE, violet, external thread, nominal volume: 4.5 ml, Cryo Performance Tested, 25 piece(s)/bag |

**Product characteristics**

|                                       |                     |
|---------------------------------------|---------------------|
| Label/ Print                          | with print          |
| Colour of print/label                 | white               |
| Graduation                            | yes                 |
| Cap                                   | QuickSeal screw cap |
| Closure type                          | screw cap           |
| Internally or externally threaded cap | external thread     |
| Volume of work                        | 4,5 ml              |
| Vessel type                           | Tube                |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Size**

|                      |          |
|----------------------|----------|
| Sample volume        | 4.5 ml   |
| Max. volume          | 5 ml     |
| Diameter             | 12.65 mm |
| Length of product    | 90.7 mm  |
| Length excluding cap | 83.8 mm  |

**Material & colours**

|                  |                                   |
|------------------|-----------------------------------|
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap    | violet                            |
| Tube material    | Polypropylene (PP)                |

**Purity & certification**

|                  |                                |
|------------------|--------------------------------|
| Product category | In vitro diagnostic device, CE |
| CE certified     | CE                             |
| Purity standard  | Cryo Performance Tested        |
| Sterilisation    | Electron irradiation           |
| Pyrogen-free     | yes                            |
| Non-cytotoxic    | yes                            |
| Batched          | yes                            |

**Packaging**

|                               |            |
|-------------------------------|------------|
| Minimum order qty.            | 250        |
| Type of smallest subpackaging | bag        |
| Piece(s) / inner box          | 250        |
| Piece (s) / outer case        | 1000       |
| Piece(s) / pallet             | 52000      |
| Depth of box                  | 293 mm     |
| Width of box                  | 183 mm     |
| Height of box                 | 118 mm     |
| Depth of case                 | 598 mm     |
| Width of case                 | 378 mm     |
| Height of case                | 135 mm     |
| Case volume                   | 0.0305 cbm |
| Weight of product             | 0.0038 kg  |
| Weight of case                | 5.03 kg    |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**CryoPure tubes, 5 ml, QuickSeal screw cap, violet**

EAN of inner box                      4038917106146  
EAN case                                4038917080699

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

