




Product Specification
Tissue culture dish, (ØxH): 100 x 20 mm, surface:
Suspension

Page 1



Product description

Order number	83.3902.500
Product description	Tissue culture dish, (ØxH): 100 x 20 mm, material: PS, surface: Suspension, for suspension cells, coding colour: green, TC Tested, 10 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Cell type	suspension cells
Base shape	SUREGrip Design
Growth area	58.8 cm ²
Volume of work	13 ml
Vessel type	Dish
Coding colour	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tissue culture dish, (ØxH): 100 x 20 mm, surface:
Suspension

Page 2

Size

Diameter	100 mm
Height of product	20 mm

Material & colours

Product material	Polystyrene (PS)
Closure material	Polystyrene (PS)

Purity & certification

Product category	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	300
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	300
Piece(s) / pallet	4200
Depth of case	505 mm
Width of case	313 mm
Height of case	460 mm
Case volume	0.0727 cbm
Weight of product	0.0193 kg
Weight of case	7.37 kg
EAN of inner box	4038917338226
EAN case	4038917338219

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

