



Product description

Order number	83.3911.502
Product description	Cell culture flask, T-75, material: PS, surface: Suspension, for suspension cells, coding colour: green, Filter cap, TC Tested, 5 piece(s)/bag

Product characteristics

Cap	Filter cap
Cell type	suspension cells
Growth area	75 cm ²
Volume of work	21 ml - 55 ml
Vessel type	flask
Coding colour	green

Material & colours

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Colour of cap	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Cell culture flask, T-75, surface: Suspension,
Filter cap

Page 2

Purity & certification

Product category	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	100
Piece(s) / pallet	1600
Depth of case	445 mm
Width of case	402 mm
Height of case	435 mm
Case volume	0.0778 cbm
Weight of product	0.05 kg
Weight of case	6.78 kg
EAN of inner box	4038917337861
EAN case	4038917366090

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT