



Product Specification
Cell culture flask, T-75, surface: Standard,
2-position screw cap

Page 1



Product description

Order number	83.3911
Product description	Cell culture flask, T-75, material: PS, surface: Standard, for Adherent cells, coding colour: red, 2-position screw cap, TC Tested, 5 piece(s)/bag

Product characteristics

Cap	2-position screw cap
Cell type	Adherent cells
Growth area	75 cm ²
Volume of work	21 ml - 55 ml
Vessel type	flask
Coding colour	red

Material & colours

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Colour of cap	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Cell culture flask, T-75, surface: Standard,
2-position screw cap

Page 2

Purity & certification

Product category	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	100
Piece(s) / pallet	1600
Depth of case	445 mm
Width of case	402 mm
Height of case	435 mm
Case volume	0.0778 cbm
Weight of product	0.0502 kg
Weight of case	6.8 kg
EAN of inner box	4038917337656
EAN case	4038917366045

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT