



### **Product description**

|                     |                                                                                                                    |
|---------------------|--------------------------------------------------------------------------------------------------------------------|
| Order number        | 83.3945.100                                                                                                        |
| Product description | Cell strainer, pore size: 100 µm, yellow, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister |

### **Product characteristics**

|                        |                                       |
|------------------------|---------------------------------------|
| Example of use         | for producing single-cell suspensions |
| Suitable for           | 50 ml centrifuge tubes                |
| Filter                 | yes                                   |
| Surface area of filter | 4.2 cm <sup>2</sup>                   |
| Pore size              | 100                                   |

### **Size**

|                   |          |
|-------------------|----------|
| Diameter          | 34.09 mm |
| Height of product | 21 mm    |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Membrane material | nylon              |
| Frame material    | Polypropylene (PP) |
| Colour of frame   | yellow             |
| Filter material   | nylon              |

### **Purity & certification**

|                  |                                                      |
|------------------|------------------------------------------------------|
| Product category | no medical device   IVD                              |
| Purity standard  | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation    | Electron irradiation                                 |
| Pyrogen-free     | yes                                                  |
| Non-cytotoxic    | yes                                                  |
| Batched          | yes                                                  |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 50            |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 50            |
| Piece (s) / outer case        | 600           |
| Piece(s) / pallet             | 9600          |
| Depth of case                 | 520 mm        |
| Width of case                 | 385 mm        |
| Height of case                | 410 mm        |
| Case volume                   | 0.0821 cbm    |
| Weight of product             | 0.0068 kg     |
| Weight of case                | 6.78 kg       |
| EAN of inner box              | 4038917415040 |
| EAN case                      | 4038917415033 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

