



**Product Specification**  
**Cell culture tubes, (LxØ): 97 x 16 mm, skirted  
conical base, TC-treated**

**Page 1**



**Product description**

Order number	83.9923.929
Product description	Cell culture tubes, 10 ml, (LxØ): 97 x 16 mm, skirted conical base, PS, TC-treated, for Suspension and mono-layer cultures, sterile, 5 piece(s)/bag

**Product characteristics**

Graduation	no
Cap	standard screw cap
Cell type	Suspension and mono-layer cultures
Base shape	skirted conical base
Volume of work	10 ml
Vessel type	Tube
Type of surface treatment	TC-treated
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Cell culture tubes, (LxØ): 97 x 16 mm, skirted**  
**conical base, TC-treated**

**Page 2**

**Size**

Diameter	16 mm
Length of product	97 mm

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Polyethylene (PE)
Colour of cap	red

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	1000
Piece(s) / pallet	20000
Depth of case	428 mm
Width of case	438 mm
Height of case	372 mm
Case volume	0.0697 cbm
Weight of product	0.0084 kg
Weight of case	9.2 kg
EAN of inner box	4038917154611
EAN case	4038917027878

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-10

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**