



Product Specification
Cell culture tubes, (LxØ): 125 x 16 mm, round base, TC-treated

Page 1



Product description

Order number	83.9923.945
Product description	Cell culture tubes, 15 ml, (LxØ): 125 x 16 mm, round base, PS, TC-treated, for Suspension and mono-layer cultures, sterile, 5 piece(s)/bag

Product characteristics

Graduation	no
Cap	standard screw cap
Cell type	Suspension and mono-layer cultures
Base shape	round base
Volume of work	15 ml
Vessel type	Tube
Type of surface treatment	TC-treated
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Cell culture tubes, (LxØ): 125 x 16 mm, round
base, TC-treated**

Page 2

Size

Diameter	16 mm
Length of product	125 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Polyethylene (PE)
Colour of cap	red

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	1000
Piece(s) / pallet	21000
Depth of case	428 mm
Width of case	438 mm
Height of case	372 mm
Case volume	0.0697 cbm
Weight of product	0.0086 kg
Weight of case	9.4 kg
EAN of inner box	4038917154635
EAN case	4038917027892

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

