




**Product Specification**  
**lumox® dish 50, Tissue culture dish, with foil**  
**base, Ø: 50 mm, Adherent cells**

**Page 1**



**Product description**

Order number	94.6077.410
Product description	lumox® dish 50, Tissue culture dish, with foil base, Ø: 50 mm, for Adherent cells, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag

**Product characteristics**

Closure type	slip-on lid
Design	with foil base
Cell type	Adherent cells
Base shape	flat base
Growth area	20.35 cm <sup>2</sup>
Volume of work	5 ml - 10 ml
Vessel type	Dish
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**lumox® dish 50, Tissue culture dish, with foil base, Ø: 50 mm, Adherent cells**

**Page 2**

**Size**

Diameter	50 mm
Height of product	12 mm

**Material & colours**

Product material	lumox®
Frame material	Polystyrene (PS)

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	50
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	200
Piece(s) / pallet	7200
Depth of case	362 mm
Width of case	162 mm
Height of case	125 mm
Case volume	0.0073 cbm
Weight of product	0.009 kg
Weight of case	4.59 kg
EAN of inner box	4038917255233
EAN case	4038917255226

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**