



Product Specification
lumox® multiwell, Cell culture plate, with foil
base, 24 well, 20 piece(s)

Page 1



Product description

| | |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number | 94.6110.024 |
| Product description | lumox® multiwell, Cell culture plate, with foil base, 24 well, frame: black, frame: PS, for Adherent cells, flat base, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister |

Product characteristics

| | |
|-----------------|---------------------|
| Closure type | slip-on lid |
| Design | with foil base |
| Cell type | Adherent cells |
| Base shape | flat base |
| Number of wells | 24 |
| Growth area | 1.9 cm ² |
| Volume of work | 500 µl - 1.500 µl |
| Vessel type | Plate |
| Coding colour | red |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**lumox® multiwell, Cell culture plate, with foil
base, 24 well, 20 piece(s)**

Page 2

Size

| | |
|-------------------|-----------|
| Width of product | 85.48 mm |
| Height of product | 23.4 mm |
| Length of product | 127.76 mm |

Material & colours

| | |
|------------------|------------------|
| Product material | lumox® |
| Frame material | Polystyrene (PS) |
| Colour of frame | black |

Purity & certification

| | |
|------------------|------------------------------------------------------|
| Product category | no medical device IVD |
| Purity standard | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation | Electron irradiation |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|----------------------------------|---------------|
| Minimum order qty. | 20 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 20 |
| Piece(s) / pallet | 1920 |
| Depth of case | 333 mm |
| Width of case | 303 mm |
| Height of case | 150 mm |
| Case volume | 0.0151 cbm |
| Weight of product | 0.065 kg |
| Weight of case | 1.67 kg |
| EAN of inner box | 4038917255448 |
| EAN case | 4038917255431 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

