



### **Product description**

Order number	83.1840.001
Product description	Coverslips, for Adherent cells, Ø: 22 mm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 25 piece(s)/bag

### **Product characteristics**

Suitable for	6-well plates
Cell type	Adherent cells
Type of surface treatment	TC-treated
Coding colour	red

### **Size**

Diameter	22 mm
----------	-------

### **Material & colours**

Product material	Polyethylene terephthalate, Glycol-modified (PET-G)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	200
Type of smallest subpackaging	bag
Piece(s) / inner box	200
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Depth of case	310 mm
Width of case	144 mm
Height of case	335 mm
Case volume	0.015 cbm
Weight of product	0.0003 kg
Weight of case	1.8 kg
EAN of inner box	4038917134040
EAN case	4038917025386

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

