



Product description

| | |
|---------------------|---|
| Order number | 83.1850 |
| Product description | Membrane pump, for pressure filtration, for Filtropur LS, 1 piece(s)/case |

Product characteristics

| | |
|--------------|-------------------------|
| Application | for pressure filtration |
| Suitable for | Filtropur LS |

Size

| | |
|-------------------|--------|
| Width of product | 99 mm |
| Height of product | 177 mm |
| Length of product | 130 mm |

Purity & certification

| | |
|------------------|--|
| Product category | Appliance, according to machinery or equipment directive, CE |
| CE certified | CE supplier |
| Batched | no |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 1 |
| Piece(s) / pallet | 50 |
| Case volume | 0.0157 cbm |
| Weight of product | 3.159 kg |
| Weight of case | 3.51 kg |
| EAN of inner box | 4038917087940 |
| EAN case | 4038917006101 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

