



Product description

Order number	83.3990.025
Product description	Screw cap, for tissue culture flask T-25, 2-position screw cap, without filter, red, TC Tested, 1 piece(s)/bag

Product characteristics

Cap	2-position screw cap
Closure type	screw cap
Suitable for	tissue culture flask T-25
Filter	no

Material & colours

Product material	High Density Polyethylene (HD-PE)
Colour of product	red
Colour of cap	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Screw cap, for tissue culture flask T-25, without filter

Purity & certification

Product category	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	100
Piece(s) / pallet	16800
Depth of case	288 mm
Width of case	218 mm
Height of case	150 mm
Case volume	0.0094 cbm
Weight of product	0.0019 kg
Weight of case	0.57 kg
EAN of inner box	4038917359924
EAN case	4038917359917

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

