




**Product Specification**  
**Serological pipette, with tip, plugged, 10 ml,  
sterile, 25 piece(s)/bag**

**Page 1**



**Product description**

Order number	86.1254.025
Product description	Serological pipette, with tip, plugged, 10 ml, coding colour: orange, material: PS, sterile, 25 piece(s)/bag

**Product characteristics**

Colour of print/label	orange
Graduation	yes
Design	with tip, plugged
Volume of work	10 ml
Coding colour	orange

**Size**

Diameter	10 mm
Length of product	335 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Serological pipette, with tip, plugged, 10 ml,  
sterile, 25 piece(s)/bag**

**Page 2**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	500
Piece(s) / pallet	27000
Depth of case	413 mm
Width of case	250 mm
Height of case	270 mm
Case volume	0.0279 cbm
Weight of product	0.0091 kg
Weight of case	5.26 kg
EAN of inner box	4038917127318
EAN case	4038917021135

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**