



## **Product description**

Order number	95.950
Product description	Felt marker, black, alcohol-resistant, application e.g.: for labelling plastic tubes and reaction tubes, 10 piece(s)/case

## **Product characteristics**

Label/ Print	with print
Colour of print/label	black
Example of use	for labelling plastic tubes and reaction tubes

## **Size**

Diameter	13 mm
Length of product	133 mm

## **Material & colours**

Product material	Polypropylene (PP)
Colour of product	white/black

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	10
Type of smallest subpackaging	case
Piece(s) / inner box	10
Piece (s) / outer case	500
Piece(s) / pallet	45000
Depth of case	145 mm
Width of case	112 mm
Height of case	58 mm
Case volume	0.0009 cbm
Weight of product	0.0161 kg
Weight of case	9.1 kg
EAN of inner box	4038917380355
EAN case	4038917380348

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

