




**Product Specification**  
**lumox® multiwell, Cell culture plate, with foil**  
**base, 96 well, 20 piece(s)**

**Page 1**



**Product description**

Order number	94.6120.096
Product description	lumox® multiwell, Cell culture plate, with foil base, 96 well, frame: black, frame: PS, for Adherent cells, flat base, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

**Product characteristics**

Closure type	slip-on lid
Design	with foil base
Cell type	Adherent cells
Base shape	flat base
Number of wells	96
Growth area	0.34 cm <sup>2</sup>
Volume of work	25 µl - 340 µl
Vessel type	Plate
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**lumox® multiwell, Cell culture plate, with foil  
base, 96 well, 20 piece(s)**

**Page 2**

**Size**

Width of product	85.48 mm
Height of product	17.2 mm
Length of product	127.76 mm

**Material & colours**

Product material	lumox®
Frame material	Polystyrene (PS)
Colour of frame	black

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	20
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	3360
Depth of case	288 mm
Width of case	218 mm
Height of case	150 mm
Case volume	0.0094 cbm
Weight of product	0.061 kg
Weight of case	1.52 kg
EAN of inner box	4038917255479
EAN case	4038917255462

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

