



**Product Specification**  
**Petri dish, 35 x 10 mm, transparent, with  
ventilation cams**

**Page 1**



**Product description**

Order number	82.1135.500
Product description	Petri dish, (ØxH): 35 x 10 mm, material: PS, transparent, with ventilation cams, sterile, 20 piece(s)/tube bag

**Product characteristics**

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	yes

**Size**

Max. volume	10 ml
Diameter	35 mm
Height of product	10 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Petri dish, 35 x 10 mm, transparent, with  
ventilation cams**

**Page 2**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Latex-free	yes
Purity standard	sterile
Sterilisation	Gamma irradiation
Safety data sheet	no
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	tube bag
Piece(s) / inner box	20
Piece (s) / outer case	500
Piece(s) / pallet	30000
Depth of case	322 mm
Width of case	257 mm
Height of case	245 mm
Case volume	0.0203 cbm
Weight of product	0.0036 kg
Weight of case	2.33 kg
EAN of inner box	4038917347839
EAN case	4038917365772

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-26

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**