




Product Specification
**Petri dish, 150 x 20 mm, transparent, with
ventilation cams**

Page 1



Product description

Order number	82.1184.500
Product description	Petri dish, (ØxH): 150 x 20 mm, material: PS, transparent, with ventilation cams, sterile, 10 piece(s)/tube bag

Product characteristics

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	yes

Size

Max. volume	275 ml
Diameter	150 mm
Height of product	20 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Petri dish, 150 x 20 mm, transparent, with
ventilation cams**

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Latex-free	yes
Purity standard	sterile
Sterilisation	Gamma irradiation
Safety data sheet	no
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	tube bag
Piece(s) / inner box	10
Piece (s) / outer case	100
Piece(s) / pallet	1800
Depth of case	438 mm
Width of case	324 mm
Height of case	540 mm
Case volume	0.0766 cbm
Weight of product	0.0524 kg
Weight of case	6.81 kg
EAN of inner box	4038917347860
EAN case	4038917347853

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-25

This document was prepared by EDP
support and is valid without signature.



SARSTEDT