




**Product Specification**  
**Petri dish, 92 x 16 mm, transparent, with  
ventilation cams**

**Page 1**



**Product description**

Order number	82.1473.001
Product description	Petri dish, (ØxH): 92 x 16 mm, material: PS, transparent, with ventilation cams, sterile, 20 piece(s)/tube bag

**Product characteristics**

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	yes

**Size**

Max. volume	80 ml
Diameter	92 mm
Height of product	16 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Petri dish, 92 x 16 mm, transparent, with  
ventilation cams**

**Page 2**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Latex-free	yes
Purity standard	sterile
Sterilisation	Gamma irradiation
Safety data sheet	no
Batched	yes

**Packaging**

Minimum order qty.	480
Type of smallest subpackaging	tube bag
Piece(s) / inner box	480
Piece (s) / outer case	480
Piece(s) / pallet	9600
Depth of case	584 mm
Width of case	387 mm
Height of case	360 mm
Case volume	0.0814 cbm
Weight of product	0.016 kg
Weight of case	8.79 kg
EAN of inner box	4038917129329
EAN case	4038917022125

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-20

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**