



Product description

Order number	68.752
Product description	Round cuvette, max. working volume: 3.5 ml, material: PP, height: 51 mm, 1,000 piece(s)/bag

Product characteristics

Volume of work	3,5 ml
----------------	--------

Size

Max. volume	4 ml
Diameter	12 mm
Height of product	51 mm

Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Round cuvette, 3.5 ml, height: 51 mm, PP

Page 2

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	100000
Depth of case	458 mm
Width of case	408 mm
Height of case	350 mm
Case volume	0.0654 cbm
Weight of product	0.0012 kg
Weight of case	6.61 kg
EAN of inner box	4038917083959
EAN case	4038917002844

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-02

This document was prepared by EDP support and is valid without signature.



SARSTEDT