



**Product Specification**  
**DishRack, height: 370 mm, yellow, for 88 petri dishes with 92 mm Ø**

**Page 1**



**Product description**

Order number	93.1646.002
Product description	DishRack, height: 370 mm, material: PC, yellow, for 88 petri dishes with 92 mm Ø, 1 piece(s)/case

**Product characteristics**

Storage capacity	88
------------------	----

**Size**

Height of product	370 mm
-------------------	--------

**Material & colours**

Product material	Polycarbonate (PC)
Colour of product	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-07

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**DishRack, height: 370 mm, yellow, for 88 petri dishes with 92 mm Ø**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      1  
Type of smallest subpackaging                      case  
Piece(s) / inner box                      1  
Piece (s) / outer case                      1  
Piece(s) / pallet                      108  
Depth of case                      383 mm  
Width of case                      193 mm  
Height of case                      190 mm  
Case volume                      0.014 cbm  
Weight of product                      1.086 kg  
Weight of case                      1.35 kg  
EAN of inner box                      4038917132565  
EAN case                      4038917030472

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

