



Product Specification
DishRack, height: 370 mm, blue, for 88 petri dishes with 92 mm Ø

Page 1



Product description

Order number	93.1646.003
Product description	DishRack, height: 370 mm, material: PC, blue, for 88 petri dishes with 92 mm Ø, 1 piece(s)/case

Product characteristics

Storage capacity	88
------------------	----

Size

Height of product	370 mm
-------------------	--------

Material & colours

Product material	Polycarbonate (PC)
Colour of product	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**DishRack, height: 370 mm, blue, for 88 petri
dishes with 92 mm Ø**

Page 2

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	108
Depth of case	383 mm
Width of case	193 mm
Height of case	190 mm
Case volume	0.014 cbm
Weight of product	0.536 kg
Weight of case	0.8 kg
EAN of inner box	4038917132572
EAN case	4038917030489

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-02

This document was prepared by EDP support and is valid without signature.



SARSTEDT