



Product description

Order number	93.1647.001
Product description	DishRack, height: 240 mm, material: PC, red, for 52 petri dishes with 92 mm Ø, 1 piece(s)/case

Product characteristics

Storage capacity	52
------------------	----

Size

Height of product	240 mm
-------------------	--------

Material & colours

Product material	Polycarbonate (PC)
Colour of product	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest case
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 1
Piece(s) / pallet 162
Depth of case 253 mm
Width of case 193 mm
Height of case 190 mm
Case volume 0.0093 cbm
Weight of product 0.617 kg
Weight of case 0.8 kg
EAN of inner box 4038917132602
EAN case 4038917030519

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

