



Product Specification

DishRack, height: 240 mm, blue, for 52 petri dishes with 92 mm Ø

Page 1



Product description

Order number	93.1647.003
Product description	DishRack, height: 240 mm, material: PC, blue, for 52 petri dishes with 92 mm Ø, 1 piece(s)/case

Product characteristics

Storage capacity	52
------------------	----

Size

Height of product	240 mm
-------------------	--------

Material & colours

Product material	Polycarbonate (PC)
Colour of product	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**DishRack, height: 240 mm, blue, for 52 petri
dishes with 92 mm Ø**

Page 2

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	162
Depth of case	253 mm
Width of case	193 mm
Height of case	190 mm
Case volume	0.0093 cbm
Weight of product	0.617 kg
Weight of case	0.8 kg
EAN of inner box	4038917132626
EAN case	4038917030533

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-02

This document was prepared by EDP support and is valid without signature.



SARSTEDT