



### **Product description**

Order number	80.625
Product description	Swabs, in the tube, no medium, (LxØ): 108 x 16 mm, stem material: PS, stem length: 83 mm, swab material: viscose, sterile, 500 piece(s)/bag

### **Product characteristics**

Label/ Print	with paper label
Closure type	screw cap
Design	in the tube
Type of medium	no medium

### **Size**

Diameter	16 mm
Length of product	108 mm
Length including cap	108 mm
Length excluding cap	101 mm
Length including cap and plunger	97 mm
Length of stem	83 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polyethylene (PE)
Colour of cap	transparent
Swab material	viscose
Stem material	Polystyrene (PS)

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	22000
Depth of case	568 mm
Width of case	348 mm
Height of case	155 mm
Case volume	0.0306 cbm
Weight of product	0.0064 kg
Weight of case	3.73 kg
EAN of inner box	4038917083003
EAN case	4038917001984

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

