



## Product Specification

Plate rotor, for 2 tubes up to 35 mm Ø and 6 tubes, up to 20 mm Ø and 6 tubes, up to 12.5 mm Ø, for SARMIX® M 2000

Page 1



### Product description

Order number	92.180.611
Product description	Plate rotor, for 2 tubes up to 35 mm Ø and 6 tubes, up to 20 mm Ø and 6 tubes, up to 12.5 mm Ø, for SARMIX® M 2000, 1 piece(s)/case

### Product characteristics

Suitable for	SARMIX® M 2000
Storage capacity	2

### Purity & certification

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-16

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Plate rotor, for 2 tubes up to 35 mm Ø and 6 tubes, up to 20 mm Ø and 6 tubes, up to 12.5 mm Ø, for SARMIX® M 2000**

**Page 2**

### Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	10
Piece(s) / pallet	10
Case volume	0.0122 cbm
Weight of product	0.2009 kg
Weight of case	2.21 kg
EAN of inner box	4038917101998
EAN case	4038917075756

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-16

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**