



**Product description**

Order number	92.180.612
Product description	Block rotor, for 40 tubes up to 11.5 mm diameter, for SARMIX® M 2000, 1 piece(s)/case

**Product characteristics**

Suitable for	SARMIX® M 2000
Storage capacity	40

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	10
Piece(s) / pallet	10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Block rotor, for 40 tubes up to 11.5 mm diameter,**  
**for SARMIX® M 2000**

Depth of case	350 mm
Width of case	240 mm
Height of case	150 mm
Case volume	0.0126 cbm
Weight of product	0.1509 kg
Weight of case	1.71 kg
EAN of inner box	4038917102001
EAN case	4038917075763

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

