



### Product description

Order number	92.181.612
Product description	Rotor, for 6 collection tubes, for tubes Ø 28 mm

### Product characteristics

Suitable for	tubes Ø 28 mm
Storage capacity	6

### Purity & certification

Product category	no medical device   IVD
Batched	yes

### Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	40
Case volume	0.0053 cbm
Weight of product	1.3595 kg
Weight of case	1.36 kg
EAN of inner box	4038917156271

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Rotor, for 6 collection tubes, for tubes Ø 28 mm**

**Page 2**

EAN case

4038917156264

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**