




Product Specification

Closure mat, for Deep Well MegaBlock® 0.5 ml
(82.1969.002) & 1.2 ml (82.1970.002 &
82.1971.002), (LxW): 79 x 121 mm, TPE,
transparent

Page 1



Product description

Order number	95.1990.002
Product description	Closure mat, for Deep Well MegaBlock® 0.5 ml (82.1969.002) & 1.2 ml (82.1970.002 & 82.1971.002), (LxW): 79 x 121 mm, material: TPE, transparent, pierceable, Free from DNA/DNase/RNase, free from pyrogens/endotoxins, 10 piece(s)/bag

Product characteristics

Suitable for	Deep Well MegaBlock® 0.5 ml (82.1969.002) & 1.2 ml (82.1970.002 & 82.1971.002)
--------------	---

Size

Width of product	121 mm
Length of product	79 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Closure mat, for Deep Well MegaBlock® 0.5 ml
(82.1969.002) & 1.2 ml (82.1970.002 &
82.1971.002), (LxW): 79 x 121 mm, TPE,
transparent**

Page 2

Material & colours

Product material	Thermoplastic elastomers (TPE)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	Free from DNA/DNase/RNase, free from pyrogens/endotoxins
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	5500
Depth of case	380 mm
Width of case	280 mm
Height of case	200 mm
Case volume	0.0213 cbm
Weight of product	0.0144 kg
Weight of case	4.1 kg
EAN of inner box	4038917135023
EAN case	4038917042048

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-25

This document was prepared by EDP
support and is valid without signature.



SARSTEDT