



Product description

Order number	86.1199
Product description	Disposal bags, volume: 24 l, (LxW): 780 x 400 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, without print, 50 piece(s)/bag

Product characteristics

Label/ Print	without print
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

Size

Width of product	400 mm
Length of product	780 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Latex-free	yes
Safety data sheet	no
Batched	no

Packaging

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	16500
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0.0182 cbm
Weight of product	0.0271 kg
Weight of case	7.05 kg
EAN of inner box	4038917085731
EAN case	4038917004367

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

