



Product description

Order number	86.1204
Product description	Disposal bags, volume: 80 l, (LxW): 1,120 x 700 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, without print, 50 piece(s)/bag

Product characteristics

Label/ Print	without print
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

Size

Width of product	700 mm
Length of product	1.120.00 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Latex-free	yes
Safety data sheet	no
Batched	no

Packaging

Minimum order qty.	150
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	150
Piece(s) / pallet	3600
Depth of case	378 mm
Width of case	122 mm
Height of case	810 mm
Case volume	0.0374 cbm
Weight of product	0.0715 kg
Weight of case	11.3 kg
EAN of inner box	4038917085786
EAN case	4038917004411

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

