



### Product description

Order number	82.1582.200
Product description	ELISA plate, 96 well, round base, without lid, material: PS, transparent, High Binding, free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic, 25 piece(s)/bag

### Product characteristics

Label/ Print	with print
Base shape	round base
Number of wells	96
Vessel type	Plate

### Size

Max. volume	310 µl
Width of product	85.2 mm
Height of product	14.35 mm
Length of product	127.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**ELISA plate, 96 well, round base, PS, transparent,**  
**High Binding**

**Page 2**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Purity standard	free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	50
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	50
Piece(s) / pallet	2000
Depth of case	279 mm
Width of case	232 mm
Height of case	365 mm
Case volume	0.0236 cbm
Weight of product	0.0344 kg
Weight of case	2.54 kg
EAN of inner box	4038917356220
EAN case	4038917356213

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**