




Product Specification
Film, free from DNase/RNase, transparent,
pierceable

Page 1



Product description

Order number	95.2092
Product description	Film, free from DNase/RNase, transparent, suitable for qPCR, pierceable, application area: -80 - 120 °C

Product characteristics

Closure type	film
Light protection	no
Suitable for qPCR	yes
Minimum temperature	-80 °C
Maximum temperature	120 °C

Size

Width of product	82 mm
Length of product	125 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Film, free from DNase/RNase, transparent,
pierceable**

Page 2

Material & colours

Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard free from DNase/RNase
Batched yes

Packaging

Minimum order qty. 100
Piece(s) / inner box 100
Piece (s) / outer case 1200
Piece(s) / pallet 57600
Case volume 0.0142 cbm
Weight of product 0.0031 kg
Weight of case 4.36 kg
EAN of inner box 4038917335676
EAN case 4038917335669

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT