



## Product description

Order number	84.1556.010
Product description	Adhesive label, (LxW): 40 x 26 mm, material: paper, white, label layout: Name + date, print: white/black, 5,000 piece(s)/roll

## Product characteristics

Colour of print/label	white/black
-----------------------	-------------

## Size

Width of product	26 mm
Length of product	40 mm

## Material & colours

Product material	paper
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	roll
Piece(s) / inner box	5000
Piece (s) / outer case	30000
Piece(s) / pallet	2880000
Depth of case	333 mm
Width of case	303 mm
Height of case	150 mm
Case volume	0.0151 cbm
Weight of product	0.0001 kg
Weight of case	2.85 kg
EAN of inner box	4038917130578
EAN case	4038917023139

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

