



**Product description**

Order number	84.2620.001
Product description	Adhesive label, (LxW): 33 x 12 mm, material: paper, white, 4,000 piece(s)/roll

**Size**

Width of product	12 mm
Length of product	33 mm

**Material & colours**

Product material	paper
Colour of product	white

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	4000
Type of smallest subpackaging	roll
Piece(s) / inner box	4000
Piece (s) / outer case	8000
Piece(s) / pallet	2592000
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0001 kg
Weight of case	1.16 kg
EAN of inner box	4038917416948
EAN case	4038917416931

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

