



**Product Specification**  
**Drop collector, length: 112 mm, PP, tube 75 x 13 mm with funnel cup**

**Page 1**



**Product description**

Order number	86.290
Product description	Drop collector, length: 112 mm, material: PP, tube 75 x 13 mm with funnel cup, funnel: HD-PE, for extracting drops for salmonella diagnosis, 250 piece(s)/bag

**Product characteristics**

Design	tube 75 x 13 mm with funnel cup
Application	for extracting drops for salmonella diagnosis

**Size**

Diameter of opening	40 mm
Length of product	112 mm

**Material & colours**

Product material	Polypropylene (PP)
------------------	--------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-03-28

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Purity & certification**

Product category	no medical device   IVD
Latex-free	yes
Safety data sheet	no
Batched	no

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	500
Piece(s) / pallet	12000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0095 kg
Weight of case	5.45 kg
EAN of inner box	4038917081085
EAN case	4038917035620

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

