



Product Specification
Stomacher bag, (LxW): 300 x 170 mm, LD-PE,
transparent

Page 1



Product description

Order number	86.9924.531
Product description	Stomacher bag, (LxW): 300 x 170 mm, material: LD-PE, transparent, Food testing laboratories

Product characteristics

Application	Food testing laboratories
Example of use	Prolonged massage in a Stomacher® machine

Size

Width of product	170 mm
Length of product	300 mm

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Stomacher bag, (LxW): 300 x 170 mm, LD-PE,
transparent**

Page 2

Purity & certification

Latex-free	yes
Pyrogen-free	no
Safety data sheet	no

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	80000
Depth of case	458 mm
Width of case	238 mm
Height of case	373 mm
Case volume	0.0408 cbm
Weight of case	9.3 kg
EAN of inner box	4038917155380
EAN case	4038917398442

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT