



Product description

Order number 91.2620.020
Product description Thermo transfer foil, black, 1 piece(s)/roll

Material & colours

Colour of product black

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 2
Type of smallest roll
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 2
Piece(s) / pallet 648
Depth of case 260 mm
Width of case 183 mm
Height of case 85 mm
Case volume 0.004 cbm
Weight of product 0.209 kg
Weight of case 0.52 kg
EAN of inner box 4038917899703
EAN case 4038917899741

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

