



Product description

| | |
|---------------------|---|
| Order number | 95.953 |
| Product description | Felt marker, blue, waterproof, application e.g.: for labelling plastic tubes and reaction tubes, 10 piece(s)/case |

Product characteristics

| | |
|-----------------------|--|
| Colour of print/label | white |
| Example of use | for labelling plastic tubes and reaction tubes |

Size

| | |
|-------------------|--------|
| Diameter | 10 mm |
| Length of product | 142 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | black |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 10 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 10 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 48000 |
| Depth of case | 300 mm |
| Width of case | 291 mm |
| Height of case | 392 mm |
| Case volume | 0.0342 cbm |
| Weight of product | 0.009 kg |
| Weight of case | 10 kg |
| EAN of inner box | 4038917084628 |
| EAN case | 4038917003391 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

