



Product description

Order number	72.1979.102
Product description	PCR plate half skirt, 96 well, transparent, High Profile, 200 µl, PCR Performance Tested, material: PP, flat deck, 5 piece(s)/minigrip bag

Product characteristics

Suitable for qPCR	yes
Base shape	conical
Construction of the plate edges	half skirt
High or low profile	High Profile
Number of wells	96
Volume of work	10 µl - 0,2 ml
Centrifugation max (RCF)	3000 U/min
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	0.2 ml
Max. volume	200 µl
Diameter	4.95 mm
Width of product	83.6 mm
Height of product	19.33 mm
Length of product	123.6 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	25
Piece (s) / outer case	100
Piece(s) / pallet	6000
Depth of case	398 mm
Width of case	358 mm
Height of case	175 mm
Case volume	0.0249 cbm
Weight of product	0.0179 kg
Weight of case	2.63 kg
EAN of inner box	4038917307253
EAN case	4038917307246

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

