



Product description

| | |
|---------------------|--|
| Order number | 72.1979 |
| Product description | PCR plate half skirt, 96 well, transparent, High Profile, 200 µl, PCR Performance Tested, material: PP, 10 piece(s)/minigrip bag |

Product characteristics

| | |
|---------------------------------|----------------|
| Suitable for qPCR | yes |
| Base shape | conical |
| Construction of the plate edges | half skirt |
| High or low profile | High Profile |
| Number of wells | 96 |
| Volume of work | 10 µl - 0,2 ml |
| Centrifugation max (RCF) | 3000 x g |
| Centrifugation temperature | 22 |
| Centrifugation time | 5 Min |
| Autoclavable (max. °C) | 121 °C |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|----------|
| Sample volume | 0.2 ml |
| Max. volume | 200 µl |
| Diameter | 4.95 mm |
| Width of product | 86 mm |
| Height of product | 19.33 mm |
| Length of product | 126.1 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 100 |
| Type of smallest subpackaging | minigrip bag |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 100 |
| Piece(s) / pallet | 5600 |
| Depth of box | 270 mm |
| Width of box | 230 mm |
| Height of box | 131 mm |
| Depth of case | 282 mm |
| Width of case | 280 mm |
| Height of case | 247 mm |
| Case volume | 0.0195 cbm |
| Weight of product | 0.0167 kg |
| Weight of case | 2.29 kg |
| EAN of inner box | 4038917570817 |
| EAN case | 4038917570800 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

