



**Product Specification**  
**PCR plate half skirt, 96 well, white, Low Profile,  
100 µl, PCR Performance Tested, PP**

**Page 1**



**Product description**

Order number	72.1982.202
Product description	PCR plate half skirt, for LightCycler® 480, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, material: PP, optimised for qPCR, 25 piece(s)/minigrip bag

**Product characteristics**

Closure type	film
Suitable for	LightCycler® 480
Suitable for qPCR	yes
Base shape	conical
Construction of the plate edges	half skirt
High or low profile	Low Profile
Number of wells	96
Volume of work	5 µl - 0,1 ml
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### Size

Sample volume	0.1 ml
Max. volume	100 µl
Diameter	5.38 mm
Width of product	85.3 mm
Height of product	15.55 mm
Length of product	127.8 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	white

### Purity & certification

Product category	no medical device   IVD
Purity standard	PCR Performance Tested
Batched	yes

### Packaging

Minimum order qty.	100
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	25
Piece (s) / outer case	100
Piece(s) / pallet	7200
Depth of case	278 mm
Width of case	208 mm
Height of case	335 mm
Case volume	0.0194 cbm
Weight of product	0.017 kg
Weight of case	2.3 kg
EAN of inner box	4038917282017
EAN case	4038917282000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

