



**Product Specification**  
**PCR single tube, 0.2 ml, PCR Performance**  
**Tested, colour mix, PP, flat cap**

**Page 1**



**Product description**

Order number	72.737.992
Product description	PCR single tube, 0.2 ml, PCR Performance Tested, red/green/yellow/blue/violet/orange, material: PP, with anti-contamination shield, flat cap, suitable for qPCR, per colour 500 piece(s)/minigrip bag

**Product characteristics**

Graduation	no
Cap	flat cap
Closure type	push cap
Suitable for qPCR	yes
Base shape	conical
High or low profile	High Profile
Number of wells	1
Volume of work	10 µl - 0,2 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**PCR single tube, 0.2 ml, PCR Performance**  
**Tested, colour mix, PP, flat cap**

**Page 2**

**Size**

Sample volume	0.2 ml
Max. volume	200 µl
Diameter	7.5 mm
Width of product	7.5 mm
Height of product	20.65 mm
Length of product	21.7 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	colour mix

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

**Packaging**

Minimum order qty.	3000
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	3000
Piece (s) / outer case	3000
Piece(s) / pallet	432000
Depth of case	323 mm
Width of case	213 mm
Height of case	145 mm
Case volume	0.01 cbm
Weight of product	0.0002 kg
Weight of case	0.76 kg
EAN of inner box	4038917120197
EAN case	4038917067676

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-25

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**