



### **Product description**

Order number	72.982.002
Product description	PCR strip of 8, up to 100 µl, Low Profile, PCR Performance Tested, transparent, material: PP, flat cap, cap enclosed, suitable for qPCR, 125 piece(s)/minigrip bag

### **Product characteristics**

Graduation	no
Cap	flat cap
Closure type	push cap
Suitable for qPCR	yes
Base shape	conical
High or low profile	Low Profile
Number of wells	8
Volume of work	5 µl - 0,1 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Sample volume	0.1 ml
Max. volume	100 µl
Width of product	16.05 mm
Length of product	76.2 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	PCR Performance Tested
Batched	yes

### **Packaging**

Minimum order qty.	125
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	125
Piece (s) / outer case	1250
Piece(s) / pallet	56250
Depth of case	300 mm
Width of case	198 mm
Height of case	495 mm
Case volume	0.0294 cbm
Weight of product	0.0016 kg
Weight of case	3.18 kg
EAN of inner box	4038917395465
EAN case	4038917395458

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

