



Product description

| | |
|---------------------|---|
| Order number | 72.985 |
| Product description | PCR plate without skirt, 12 strips of 8 tubes (72.985.002) preloaded in PCR work tray (95.987.002), 96 well, transparent, High Profile, 200 µl, Biosphere® plus, material: PP, 1 piece(s)/blister |

Product characteristics

| | |
|---------------------------------|---|
| Closure type | push cap |
| Design | 12 strips of 8 tubes (72.985.002) preloaded in PCR work tray (95.987.002) |
| Suitable for qPCR | yes |
| Base shape | conical |
| Construction of the plate edges | without skirt |
| High or low profile | High Profile |
| Number of wells | 96 |
| Volume of work | 10 µl - 0,2 ml |
| Centrifugation max (RCF) | 3000 x g |
| Centrifugation temperature | 22 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

PCR plate without skirt, 96 well, transparent, High Profile, 200 µl, Biosphere® plus, PP

Page 2

Centrifugation time 5 Min
Autoclavable (max. °C) 121 °C

Size

Sample volume 0.2 ml
Max. volume 200 µl
Diameter 5.4 mm
Width of product 79 mm
Height of product 20.7 mm
Length of product 115 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard Biosphere® plus
Sterilisation Gas sterilisation
Pyrogen-free yes
Batched yes

Packaging

Minimum order qty. 20
Type of smallest subpackaging blister
Piece(s) / inner box 20
Piece (s) / outer case 20
Piece(s) / pallet 3360
Depth of case 288 mm
Width of case 218 mm
Height of case 150 mm
Case volume 0.0094 cbm
Weight of product 0.0228 kg
Weight of case 0.76 kg
EAN of inner box 4038917084840

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**PCR plate without skirt, 96 well, transparent, High
Profile, 200 µl, Biosphere® plus, PP**

Page 3

EAN case

4038917003612

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT