



Product description

Order number	72.994.992
Product description	PCR strip of 8, High Profile, PCR Performance Tested, red/green/blue/violet/yellow, material: PP, round cap, cap attached

Product characteristics

Graduation	no
Cap	round cap
Closure type	push cap
Base shape	conical
High or low profile	High Profile
Number of wells	8

Size

Sample volume	0.2 ml
---------------	--------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	colour mix

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	1200
Piece(s) / inner box	120
Piece (s) / outer case	1200
Piece(s) / pallet	1200
Case volume	0.0299 cbm
Weight of product	0.0031 kg
Weight of case	4.3 kg
EAN of inner box	4038917124331
EAN case	4038917080743

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

