



Product description

Order number	95.1992
Product description	PCR film, for PCR plates (REF 72.1978/.010, 72.1979/.010/.003/.201/.700, 72.1980/.010/.201/.600/.700, 72.1981/.010) & PCR strips (REF 72.985.002/.092/.992), free from DNase/RNase, transparent, suitable for qPCR, self-adhesive, application area: -80 - 120 °C, 100 piece(s)/box

Product characteristics

Closure type	film
Light protection	no
Suitable for	PCR plates (REF 72.1978/.010, 72.1979/.010/.003/.201/.700, 72.1980/.010/.201/.600/.700, 72.1981/.010) & PCR strips (REF 72.985.002/.092/.992)
Suitable for qPCR	yes
Minimum temperature	-80 °C
Maximum temperature	120 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard free from DNase/RNase
Batched yes

Packaging

Minimum order qty. 100
Type of smallest box
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 5000
Piece(s) / pallet 90000
Depth of box 90 mm
Width of box 46 mm
Height of box 170 mm
Depth of case 380 mm
Width of case 380 mm
Height of case 380 mm
Case volume 0.0549 cbm
Weight of product 0.0036 kg
Weight of case 20 kg
EAN of inner box 4038917940078
EAN case 4038917940085

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

