



## Product description

Order number	95.1994
Product description	PCR film, free from DNase/RNase, material: PO, transparent, suitable for qPCR, self-adhesive, application area: -70 - 110 °C, strong adhesion, high chemical resistance, 100 piece(s)/bag

## Product characteristics

Closure type	film
Light protection	no
Suitable for qPCR	yes
Minimum temperature	-70 °C
Maximum temperature	110 °C

## Size

Width of product	79 mm
Length of product	141 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polyolefin (PO)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Purity standard	free from DNase/RNase
Batched	yes

**Packaging**

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1200
Piece(s) / pallet	288000
Case volume	0.0052 cbm
Weight of product	0.002 kg
Weight of case	2.72 kg
EAN of inner box	4038917088114
EAN case	4038917044301

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

